

New firmware version 01.06.02 for your CoaguChek INRange / CoaguChek Vantus (US only) meter

Date:

Nov 2020

Rev. 1.0

We herewith introduce the new firmware version 01.06.02 for your CoaguChek INRange / CoaguChek Vantus (US only) meter.

This software version cannot be downgraded.

Please note:

The difference between the former firmware version 01.05.00 and the new firmware version 01.06.02 is explained as follows:

- Improvement of stability of Bluetooth LE interface

Should you have any question about the CoaguChek INRange / CoaguChek Vantus (US only) firmware version 01.06.02, please contact your Roche Service organization.

Sincerely